

When taking your impressions, please be sure to capture the following anatomical landmarks:

Edentulous Maxillary Cast

- Rugae
- Vestibule
- Tuberosity
- Hamular Notch

Edentulous Mandibular Cast

- Retromolar Pad
- External Ridge
- Labial | Lingual Sulcus


Dentulous Maxillary Cast

- Hamular Notch
- Gingival Roll
- Vestibule

Trimming Requirements

- 3mm Ledge
- Hamular Notch
- 3mm Vestibule
- 15mm Base

Maxilla



Mandible

