

## Full Denture

- Photos (see photo series)
  - Alginate impression to make custom tray
  - Final impression using custom tray
  - Wax-rims
  - Tooth shade
- 

## Nightguard

- Upper and lower impression
  - Bite registration
  - Type of nightguard and which arch
- 

## Repair

- For partial denture, pick up impression with partial in place with opposing model
  - For break, just the appliance is needed
  - If adding a new tooth, need tooth shade
- 

## Flipper

- Impression
  - Bite
  - Tooth shade
  - Are existing root tips being extracted?
- 

## Clear Duplicate

- Patient's current denture

## Partial Denture

- Photos (see photo series)
  - Alginate impression to make custom tray
  - Final impression using custom tray
  - Tooth shade
  - Cast frame or Duraflex (Flexible)
- 

## Reline

- Reline impression  
(2mm thick, relieve denture base if necessary)
- 

## Custom Tray

- Preliminary alginate impression capturing all borders and necessary landmarks including the hamular notch and retro molar pad
- 

## Surgical Guide

- Dicom, final impression or intraoral scan
- 

## Implant Bar

- Different types of bar designs have different initial requirements. Please contact us for specifics at 480-446-7063.