

Full Denture

- Photos (see photo series)
 - Alginate impression to make custom tray
 - Final impression using custom tray
 - Wax-rims + bite
 - Tooth shade
-

Immediate Denture

- Photos (see photo series)
 - Final Impression
 - Bite Registration
-

Nightguard

- Upper and lower impression
 - Bite registration
 - Type of nightguard and which arch
-

Repair

- For partial denture, pick up impression with partial in place with opposing model
 - For break, just the appliance is needed
 - If adding a new tooth, need tooth shade
-

Flipper

- Impression
 - Bite
 - Tooth shade
 - Are existing root tips being extracted?
-

Partial Denture

- Photos (see photo series)
 - Alginate impression to make custom tray
 - Final impression using custom tray
 - Opposing impression + bite
 - Tooth shade
 - Cast frame or Duraflex (Flexible)
-

Reline

- Reline impression
(2mm thick, relieve denture base if necessary)
-

Custom Tray

- Preliminary alginate impression capturing all borders and necessary landmarks including the hamular notch and retro molar pad
-

Surgical Guide

- Dicom, final impression or intraoral scan
-

Implant Bar

- Different types of bar designs have different initial requirements. Please contact us for specifics at 480-446-7063.
-

Clear Duplicate

- Patient's current denture
-