

When taking your impressions, please be sure to capture the following anatomical landmarks:

## Edentulous Maxillary Cast

- Rugae
- Vestibule
- Tuberosity
- Hamular Notch

## Edentulous Mandibular Cast

- Retromolar Pad
- External Ridge
- Labial | Lingual Sulcus


## Dentulous Maxillary Cast

- Hamular Notch
- Gingival Roll
- Vestibule

## Trimming Requirements

- 3mm Ledge
- Hamular Notch
- 3mm Vestibule
- 15mm Base

## Maxilla



## Mandible

