



Implant Torque Specification Guide



WIAND DENTAL
LABORATORY

Ensure you are applying the right amount of torque for your patient's implant restoration.

Implant	Torque (Ncm)	Torque Settings	Driver & Head
BIOHORIZONS			
Internal 3.5, 4.5, 5.7	30 Ncm	0.050" (1.25 mm) Hex	
Tapered Internal 3.5, 4.5, 5.7	30 Ncm	0.050" (1.25 mm) Hex	
External Implant 3.5	30 Ncm	0.050" (1.25 mm) Hex	
External Implant 4.0, 5.0	30 Ncm	0.050" (1.25 mm) Hex	
Straight Multi-Unit Abutment	30 Ncm	0.050" (1.25 mm) Hex	
Angled Multi-Unit Abutment	30 Ncm	0.050" (1.25 mm) Hex	
Prosthetic Screw	15 Ncm	0.050" (1.25 mm) Hex	
BIOMET 3I			
Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	20 Ncm	0.048" Hex	
Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	20 Ncm	0.048" Hex	
Certain™ XP 4/5, 5/6	20 Ncm	0.048" Hex	
Miniplant™ 3.25	20 Ncm	0.048" Hex	
Osseotite™ 3.75, 4.0, 5.0, 6.0	20 Ncm	0.048" Hex	
Osseotite™ XP 4/5, 5/6	20 Ncm	0.048" Hex	
Straight Low Profile Abutment	20 Ncm	0.048" Hex	
Angled Low Profile Abutment	20 Ncm	0.048" Hex	
Prosthetic Screw	10 Ncm	0.048" Hex	
DENTSPLY			
ANKYLOS® C/X - C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	15 Ncm	1.0 mm Hex	
ANKYLOS® C/X - /X 3.5, 4.5, 5.5, 7.0 (indexed)	15 Ncm	1.0 mm Hex	
OsseoSpeed™ EV 3.0 - Green	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 3.6 - Purple	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 4.2 - Yellow	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 4.8 - Blue	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 5.4 - Brown	25 Ncm	0.050" Hex	