

# Digital Scanning for Assistants

How to Accurately & Adequately  
Capture Data for Dental Restorations



WIAND DENTAL  
LABORATORY

**COURSE DESCRIPTION:** Starting with an accurate and adequate scan will help you and your doctor to set up every case for success *while the patient is in your chair!* Now the standard of care, intraoral scanning is finally viable for all types of restorations. Digital scans benefit all parties of the dental team as well as the patients we treat together in so many ways! It is critical to know that a digital record requires you to capture the same data as a conventional impression which is different for C&B versus full-arch. With C&B you are simply scanning for tooth structure. However, full arch requires you to capture much more data including anatomical landmarks. Join us for a thorough review, demonstration and **HANDS-ON** practice on how to capture all the data needed for each type of dental restoration and become the scanning expert at your practice!

## LEARNING OBJECTIVES:

- Explore the benefits of digital scanning
- Gain an understanding of why different data is required for each type of restoration
- Review landmarks necessary for full arch scans
- Review scanning protocols from start to finish for both C&B and full arch
- Discover scanning techniques for restorative success

Many thanks to iTero for  
sponsoring our event!



## EVENT DETAILS

**Date:** Friday, October 18th, 2024  
**Time:** 8:30am – 12:30pm (8:00am check-in/breakfast provided)  
**Location:** Wiand Dental Lab  
10555 N. 114<sup>th</sup> Street – Suite 101 – Scottsdale  
**CE:** 4 credit hours      **Fee:** \$98 to be applied as lab credit

SCAN QR CODE TO REGISTER!



SEATING IS  
LIMITED RESERVE  
YOURS TODAY!

### PRESENTER:

Scott Lowrey  
Technical Education  
and Support Manager  
Wiand Dental Lab, Inc.

